Extrication/spine immobilization device Extricator/inmovilizador espinal

This appliance conforms with the Directive 93/42/CEE "Medical Devices".

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All immobilizing procedures must be carried as prescribed by the local protocols in use and following the local specifications. Immobilize the patient's head with arms, in order to maintain the spinal column alignment.

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Apply a rigid cervical collar of the suitable size.

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Position the device behind the patient's back and align it to the spinal column with the belts facing towards the outside.

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Once the SED is in the correct position, before fixing it around the patient's body, pull leg straps from the rear and leave them loose on sides. Wrap chest supports around patient's chest. Adjust device position by lifting it until the abdominal flaps support armpits.

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Buckle the thorax belts according to the following sequence generally considered correct: • Intermediate (yellow) • Lower (red) • Upper (green)
This sequence can not be indicated in all rescue situations. Follow your local EMS guidelines. If the device is used for extrication of pregnant women, it is necessary to fold the thoracic flaps offering more space to the abdomen.
One possible way to close the device, variable depending on the breast and abdomen size, is the following:
• Intermediate belt in the corresponding buckle • Lower belt in the upper buckle
• Upper belt in the lower buckle

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Close the leg straps sliding them under the patient's legs and attach them to the appropriate buckles located on the opposite side.

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Tighten both leg straps during the patient's transfer. In case of suspected femoral or pelvic fractures, attention must be paid when using the thigh belts.
The handles of the device can help to achieve better vertical positioning of the device and its best alignment.

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Check that the SED is correctly fitted according to patients's size and proceed to immobilize head. Don't cause neck and head movements that could compromise the patient's clinical status. If necessary use one of special cushions to fill in any gap between neck and device. Fix the chin strap making sure that the cervical collar is supported and inclining it towards the patient's ears.

This will avoid having to immobilize the jaw which means the patient is unable to open his mouth. Fix the forehead strap perpendicular to the surface of the forehead bone. All straps must be fixed symmetrically.

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Once the correct closure of all fixings and the condition of the patient have been checked, the patient can be moved. To move the patient the thorax belt must be tightened. It should then be loosened when the patient is fixed to the spinal board. During use, pressure and lengthwise flexion must not exceed 12 kg. If the patient is vertically lifted, make sure that the force is equally distributed on all handles according to the maximum loading capacity of the device. When the patient has been transferred to another device such as spinal board or a stretcher, immobilize him using proper fixing devices.

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• The use of the device requires the presence of at least two trained and qualified operators.
• Avoid torsions of the patient's thorax during the application of the device and belts adjustment.
• The use with other devices or systems that are not approved by the manufacturer, may cause serious injury to the patient and could cause accidents.

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SR00111B SED Spencer
Extrication Device

SR00102B SED XS

SR00100B SED XL

User's Manual

1 GENERAL INFORMATION
1.1 Aim and contents <p>The aim of this manual is to supply all the information necessary so that the client, will not only attain adequate use of the appliance, he will also be capable of using the instrument in the most autonomous and secure way possible. This includes information regarding technical aspects, functioning, maintenance, spare parts and safety.</p> 1.2 Conservation of the instruction manual <p>The instruction and maintenance manual must be kept with the product, inside the specially provided container and above all, away from any substances or liquids which could compromise perfect legibility.</p> 1.3 Symbols used
 General or specific warning
 See instructions for use
 Lot number
 Product code
 Product compliant with specifications of the Directive 93/42/CEE

1.4 Servicing request

For any information regarding the use, maintenance and installation, please contact the Spencer Customer Care Service on +39 0521 541111 - Fax +39 0521 541222 e-mail: export@spencer.it for general informations about Spencer products, service@spencer.it for technical assistance or write to Spencer Italia S.r.l. - Via Provinciale, 12 43038 Sala Baganza (Parma) ITALY.

1.5 Demolition

Follow the current regulations.

1.6 Labelling

Each device has got an identifying label, positioned on the device itself and/or on the box. This label includes information about the Manufacturer, the product, CE mark, lot number (LOT) or serial number (SN).

2 WARNINGS

2.1 General warnings

- Before carrying out any kind of operation on the appliance (training, installation, use), the operator must carefully read the enclosed instructions, paying particular attention to the correct safety precautions and to the procedures to be followed for installation and for correct use.
- In the case of any doubts as to the correct interpretation of the instructions, please contact Spencer Italia S.r.l. for any necessary clarifications.
- Regularly check the appliance.
- In the case of any abnormalities or damage to the appliance, which could jeopardize the functioning and the safety, the appliance must be immediately removed from service.
- Do not alter or modify in any way the appliance; any such interference could cause malfunctions and injury to the patient and/or rescuer.
- The appliance must not in any way be tampered with. In such cases all responsibility will be denied for any malfunctions or injuries caused by the appliance itself.
- Those who modify or have modified, prepare or have prepared medical appliances in such a way that they no longer serve the purpose for which they were intended, or no longer supply the intended service, must satisfy the valid conditions for the introduction onto the market.
- Ensure that all the necessary precautions are taken in order to avoid the hazards that can arise as the result of contact with blood or body fluids.
- Handle with care.
- Register and store with these instructions: lot number, place and date of purchase, first date of use, date of checks, name of users, any comments.
- Use of the device in anyway other than described in this manual is forbidden.

2.2 Specific warnings

- The product must be used by trained personnel only.
- If any failure or incorrect functioning of the device is detected, it must be immediately substituted with a similar item so that the rescue procedures are guaranteed without any interruption.
- When the device is being used, the assistance of qualified staff must be guaranteed.
- The extrication/spine immobilization device should not be exposed nor get in contact with heat sources or flammable agents.
- Store in a cool, dry, dark place and do not expose to direct sun. Do not store the device underneath any heavy objects.
- Use the device only as described in this manual.
- Always check the integrity of all the parts of the device before use.
- Replace immediately damaged parts.
- Do not flex the SED in excessive and extended way.
- Device must be carried out in compliance with the local health authority regulations.
- The device must be used following the procedures approved by the local Emergency Medical Services.

2.3 Contraindications and side effects

The use of this device, if used as described in this manual, does not present any contraindications or collateral effects.

3 PRODUCT DESCRIPTION

3.1 Intended use

SED Spencer Extrication Device is an instrument for extrication and spine immobilization, ideal for all emergency situations. Particularly suitable for victims blocked in vehicles, collapsed buildings and uneasy places. It can also be used with pregnant women, children and babies. It is ideal for the immobilization of the hips and pelvis fractures. The vertical rigidity and the horizontal flexibility are the main characteristics of this device, which allows a rapid immobilization of the spinal column and a secure extrication. SED is completely X-ray compatible. The patient cannot intervene on the device in any way.

3.2 Main components (fig. 11)

- A Covering made of 400 D nylon covered in vinyl
 - B Internal supports made of birch splints
 - C Head support
 - D Strips of self-adherent strips to fix the nape and chin straps
 - E Abdominal buckles with 50 mm polypropylene tape
 - F Abdominal belts with 50 mm polypropylene tape
 - G Tight buckles with 50 mm polypropylene tape
 - H Tight belts with 50 mm polypropylene tape
 - I Lifting handles
 - J Nape strap made of 30 mm tape
 - K Chin strap made of 30 mm tape
 - L Posterior nape cushion and support filled with indestructible expanded polyurethane
 - M Cross-sectional strap
 - N Transport bag
- 3.3 Models**
- See back.
- 3.4 Technical data**
- Dimensions folded inside transport bag: 880 x 110 x h305 mm
Width extended: 900 mm
Length: 830 mm
Weight with bag and accessories: 3 kg
Load capacity: 230 kg
Functioning temperature: from -20 to +50 °C
Storage temperature: from -20 to +50 °C
Relative humidity: maximum 50%

4 OPERATING INSTRUCTIONS

4.1 Transport and storage

Before transporting the appliance, make sure that it is correctly packed ensuring also that there are no risks of shocks, bumps or falls during the transport itself.

Keep the original packaging for use in case of any further transport. Damage to the appliance caused during transport and handling is not covered by the guarantee. Repairs or replacement of the damaged parts are the responsibility of the client. The device must be stored in a dry, cool area away from direct sunlight. It must not be placed in contact with any substances or chemical agents which could cause damage and reduce safety characteristics.

4.2 Preparation

On receipt of the product:

- remove the packaging and display the material so that all components are visible
- check that all the components/pieces on the accompanying list are present
- the appliance must be checked before every use so as to reveal any working abnormalities and/or damage caused by transport and/or storage. If you want to check the device, all straps on the SED must be correctly attached so as to test using all straps and not just one buckle and strap. Using just one buckle and strap could lead to abnormal strain which would not be a conform simulation of the normal use.

- Before putting the device in service make the following careful checks:
- general integrity (absence of cuts, holes, abrasions)
 - self-adherent strips and belts tightness
 - integrity of splints (apply light pressure on the centre of the device in order to check any breakages; do not apply a flexion force higher than 12 kg)

4.3 Functioning

See back.

4.4 Troubleshooting

Excessive vertical flexion

The belt buckle will not attach correctly

Cause

Probably one or more of the splints are broken

Probable buckle break

Remedy

Put immediately the device out of service and contact the service centre

5 MAINTENANCE AND CLEANING

5.1 Cleaning

The operator must always wear adequate personal protection such as gloves and mask etc. during all checking and cleaning procedures. Failure to carry out cleaning operations can cause the risk of cross infection. For a correct device storage, clean it as follows. Clean by means of a clean cloth and bactericidal or germicidal disinfectant. Rinse with hot water, extend the device and leave it to dry in a parked environment free from humidity. Keep clean the outward surface by means of a clean cloth. Do note use chemical solvents

5.2 Maintenance

Establish a maintenance programme and periodic testing, identifying a reference employee.

The inspection frequency is determined by factors such as legal requirements, the type of use, frequency of use, environmental conditions during use and storage.

The person responsible for every day maintenance can only substitute the spare parts indicated on paragraph 6.2 "Spare Parts". All other substitutions or repairs can be carried out only by the manufacturer or by a centre authorised by the manufacturer.

For any operations that are not carried out directly by the manufacturer but by an authorised centre, we have to underline that a report regarding all operations carried out must be requested. This will permit both Spencer Italia S.r.l. and the end user to keep a log book regarding the operations carried out on the device. The device, if used as indicated in the following instruction manual, has an average life span of 5 years.

5.2.2 Special servicing

Only the manufacturer or centres with written authorisation are authorised to complete any special servicing operations. For any operations that are not carried out directly by the manufacturer but by an authorised centre, we have to underline that a report regarding all operations carried out must be requested. This will permit both Spencer Italia S.r.l. and the end user to keep a log book regarding the operations carried out on the device. The device, if used as indicated in the following instruction manual, has an average life span of 5 years.

6 ACCESSORIES AND SPARE PARTS

- 6.1 Accessories**
 - SR00101B SED XL Kit
 - SR00220A Transport bag
- 6.2 Spare parts**
 - SR00210A Nape and chin straps replacement kit.
 - SR00221A Posterior cushion with self-adherent strips

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All rights reserved. No part of this document can be photocopied, reproduced or translated into another language without the written approval of Spencer Italia S.r.l.