

Since Today the Fixation System for Ambulance has been changed Forever



The New Standard For Fixation of Medical Equipments In Ambulance

The safety of an ambulance is measured by the capacity to integrate various medical equipment in perfect harmony and synergy in both simulations and under actual testing.

The Spencer fixation systems have been tested and optimized in terms of crash behaviour, robustness, rigidity as well as resistance to bending and torsion.

- * High power compact Suction unit fixation system for all types of spencer suction machines .
- * Integrated 10 G fixation mounts for different types of ventilators in a ambulance .
- * Oxygen tank holder with 20 G Fixations .



- * Fixation system for first aid boxes .
- * Fixation device for emergency kits and resuscitation kits in a ambulance .
- * 10 G and 20 G fixation system for ambulance chairs and evacuation chairs.



Oxypack

The versatile & compact fixation to attached oxygen cylinder to stretcher & trolley in a ambulance.



Fix Board

Wall mounted fixation system for all types of spine board in ambulance

The Products of Innovation



SVX / SXO
Vertical Fixation System for Scoop Stretcher

It allows you to vertically or horizontally fix the scoop stretcher on the wall of the ambulance. Made of aluminium and glass loaded nylon, it guarantees a fast and safe positioning.



F 400
20 G / 10 G Fixation system

A tested 20 G blocking system, dedicated to the Spencer chairs which permits an easy, fast and safe fixture of the chairs to the ambulance surface or to special supports.



G- Max
20 G Fixation System

A fixation system which allows you to fix all types of self loading ambulance stretchers into different types of ambulance platform.