

The world's biggest EMS catalogue

NOW IN INDIA

The leading European producer of ambulance equipment, Spencer has been inventing solutions for emergency services since 1989. Functionality, safety, long span and performance are the keys of a range of more than 1700 products.

Now this range is closer to you.

India is one of the fastest growing market for EMS, and over the past few years, India professionals have become very discerning and very demanding in terms of affordable and reliable devices.

Spencer India is the answer to this request.

Spencer in the answer to emergency.

www.spencer.it



LOCAL CONSULTANCY



LOCAL WAREHOUSE



LOCAL SUPPORT SERVICES

SPENCER OXYGEN DELIVERY SYSTEMS

Spencer offers complete range of oxygen delivery for type A, B, C and D ambulances.

Spencer newly designed and customized oxygen delivery solutions are ideal solution for modern oxygen solutions for advance ambulance.

Spencer offers all components for oxygen solution for different type of ambulance.



02 Complete Solution for Ambulance

Innovative design, functional technology and flexibility in use. The Spencer system for oxygen, compressed air and vacuum, modular and flexible, can be tailor-made so as to respond precisely to each and every specific use. A unique and complete solution for the distribution of medical gas on the ambulance.

Especially designed for National Ambulance code AIS-125.

POP Modular wall system for medical gas distribution

Complete, modular and tailor-made

Modules with different operational and diagnostic specializations

Fast and safe installation, easy to clean

Extruded compensation chamber in plug rack

Conform to UNI EN-1789:2010, UNI EN ISO 9170-1 and 10 G tested



Rek-OX Modular external wall support pipeline for medical gas

Modular and tailor-made

Compensation chamber extruded in bar

Compatible with the exchange valve Self Switch and the led display Monitox

Possibility to fix the clips Rek and Rek Pro directly on the bar

Conform to UNI EN-1789:2010, UNI EN ISO 9170-1 and 10 G tested

Self Switch Exchange valve for oxygen systems

The only one-piece exchanger on the market

Changes automatically from one tank to the other

Autonomous function without the need of electrical supply

Fast and safe installation by using the holes placed on the aluminium structure

Conform to the directive 93/42/EEC



Monitox LED display for oxygen tanks content

Visualization of the oxygen content in the tanks also when mounted in a position difficult for the operator to reach

Easy installation with Velcro®

Coded and dedicated connections for correct installation

Composed of an easily installed Kit

Conform to the European directive 93/42/EEC



**PRX
Diaphragm pressure regulator**

- Possibility to fix directly to the tank
- Innovative, efficient and secure construction
- Pressure control of the tanks can be read with a precision gauge
- Medical gas inlet connection as required by the European standards
- Conform to the directive 93/42/EEC and the norms EN 10524-1, EN 562



**FL-X
FLOWMETER/HUMIDIFIER GROUP**


- The humidifier helps to reduce the dryness of mucosas during administration of oxygen.
- Hose output for the connection of devices used for oxygen administration
- Humidifier intended to be connected to a flowmeter within an oxygen distribution systems.

MANUAL AND AUTOMATIC CHANGEOVER OXYGEN SYSTEMS FOR AMBULANCE

- We provide ambulance builders with a tailor made installation for all standards and specifications worldwide.
- The Spencer oxygen delivery system have all the components required for the oxygen delivery in a ambulance as: Pressure Regulators for each of the supply sources(stationary and portable) Low pressure, electrically conductive, Hose for oxygen.
- Oxygen piping unexposed and secured against damages, and accessible for inspection.
- 1/2/3 No. of Pressure Regulator as Required.
- The display panel has 3 values displayed to indicate the pressure of both the cylinders and distribution pressure level.
- Digital/Led Display for Patient Cabin.
- The oxygen supply comprise of pressure regulators and terminal units or flow metering devices.
- Flow metering devices for terminal units are of dial type & and without any floats.
- The Outlets have hexagonal profile to allow only matching adapters.
- The pressure regulators are directly connected to the source of supply.



VACUUM OUTLET



G-RAIL



OXYGEN INLET



MONITOX

CYL 1	0075.7	bar ²
CYL 2	0107.7	bar ²
OUTPUT	0004.4	bar ²

DIGITAL DISPLAY



TANK UP



TANKER HOLDER



OUTLET



**VACUUM
REGULATOR**



**VENTURI
ASPIRATOR**



**ODS PANEL
WITH CHANGEOVER**

Please contact our local representative for more Information about oxygen delivery system. Spencer also offer instillation & technical support for oxygen delivery system to ambulance builders.



Visit our website www.spencer-india.com

SPENCER india

Follow us



Spencer India Technologies Pvt Ltd



SPENCER INDIA TECHNOLOGIES PVT. LTD
PLOT No. 477, GROUND FLOOR, SECTOR 26, BHARTHAL,
110077 | DWARKA | NEW DELHI | INDIA

Mobile: +91-99715 86305
Email: info@spencer-india.com