




# TOTAL SAFETY SOLUTIONS



FROM COVID -19



## PRODUCT RANGE

### NITRILE EXAMINATION GLOVES



- Nitrile gloves protect from viruses and bacterial infections.
- Manufactured from synthetic nitrile.
- Unique purple color separates them from latex and synthetic gloves
- Color: Purple and Violet
- Size: Large, Medium, Small
- Powder-Free and powdered
- Gloves have very low levels of chemical allergens.
- Packing - Box 100 pcs, Box 50 Pcs / 20 Pcs / 10 Pcs / 2 Pcs

### LATEX EXAMINATION GLOVES



- Powdered Smooth Latex Examination Gloves; Color: White
- Single use disposable for multi-purpose
- Protect workers from a broad range of infections
- Unique white color distinguishes them from other types of gloves
- Size: Large, Medium, Small
- Powder-Free and powdered
- Gloves have very low levels of chemical allergens.
- Packing - Box 100 pcs, Box 50 Pcs / 20 Pcs / 10 Pcs / 2 Pcs

### STERILE LATEX GLOVES



- Sterile Latex Gloves, multipurpose for safety against infections
- Material: Powdered Smooth Latex sterile Gloves; Color: White
- Single use disposable for multi-purpose
- Sterile; package includes:
- Size: Large, Medium, Small
- Powder-Free and powdered
- Gloves have very low levels of chemical allergens.
- Packing - Box 100 pair, Box 50 pair / 20 Pair / 10 Pair / 2 Pair

### PLASTIC DISPOSABLE GLOVES



- Disposable Gloves / Transparent / Soft And Clean
- Economic In Use / Soft, Clean
- Can Be Use In Everyday use purpose
- Make yourself Safe From Bacteria.
- Can be used in Home, office, Hotels, Catering, Hospital, Pharma, Agriculture, Textile.
- Economic in Use. Soft, clean, Lint free and resistance to fungal growth
- Pack - 100 pcs, Pack - 50 Pcs, Pack - 2 pcs



## PRODUCT RANGE

### RUBBER SAFETY GLOVES



- Soft and Fluffy lining inside the Gloves, Silky surface and powder free interior
- Applications: Household chores, sanitization, cleaning & janitorial.
- Popular general purpose glove,
- Specially formulated rubber compound
- Natural rubber minimal latex protein level, Improved Fat and Abrasion Resistance
- Material: Natural Rubber
- Weight : 80 grams approx.
- Packing:- 1 Pair

### ADULT FACE SHIELD



- Designed for Adult Face – complete face covered
- Splash, water and dust resistant, disposable, easy to clean
- designed to protect your whole face from spray and splatter; droplet, dust
- General protection, suitable for home, shop or dental use
- Lightweight & Comfortable
- Equipped with an elastic band and sponge headband
- High-quality optically and distortion-free film.
- Clear for maximum visibility.
- Available in sizes - 0.3 microns, 0.4 micron, 0.5 micron, 1 micron

### CHILD / KIDS FACE SHIELD



- Specially designed for Child / Kids
- Face Shield for Children with Kids friendly Design
- Large clear visor; to maximum visibility, durable and practical.
- Equipped with an elastic band and forehead sponge, easy to adjust and install.
- Make it comfortable to wear a long time
- Front PE Sheet may have School name or Child Name
- Pack of 10 / 5/ 3 /1 Shields

### DESIGNER FACE MASK - PRINTED



- Face Mask Anti-pollution Dustproof, Washable Face Cover Men Women Unisex
- Front Cover - Poly canvas, Inner layer - Non-woven filter cloth
- TRIPLE LAYERED MASK: The Top layer is of Poly Canvas, mid layer of Non-Woven
- After every use, this mask should be soaked and washed for 15 minutes
- The mask is reusable up to 20/30 washes.
- The mask has a unique and attractive design & avoid rashes or itching.
- mask is equipped with elastic loops that allow it to fit properly and keep it in place
- The print quality is good & the mask is designed to provide you comfort

## PRODUCT RANGE

### VIRUS BLOCK OUT STERILIZATION CARD



- Sterilization Card, Home Disinfection Card with Air Sterilization Card
- exquisite appearance, wear it around your neck or put it in your pocket, The chlorine contained in the product elevator can eliminate space bacteria.
- It is suitable for, stations, schools, hospitals, offices, and kill the Viruses etc
- The product can be hung above the neck or other places where sterilization is needed.
- After the product is opened, the chlorine dioxide antibacterial component contained in the product begins to clean up and eliminate viruses and bacteria in the surrounding floating air.
- The product is valid for about 30 days, Card Size: 100\*75\*2mm, Materials: Sodium Chlorite

### CORONA SIGNAGE KIT



- Coronavirus COVID-19 Signs
- Social Distancing A-Frame Signs
- Social Distancing Banners
- Social Distancing Emoji Signs and Social Distancing Floor Shapes
- Floor Signs - Social Distancing and Social Distancing Floor Tapes
- Social Distancing Labels and Social Distancing Placards
- Posters - Social Distancing and related to Covid - 19
- Social Distancing Strips
- Social Distancing Wall Signs

### AUTOMATIC SANITIZER DISPENSER



- It is designed to nebulize hand sanitizer or alcohol, and provide automatic dosage of spray, which eliminates cross contamination, achieving optimal hand hygiene.
- automatic hand sanitizer dispenser of 2/4/5/8/10 litre ,suitable for use of disinfectant liquids, and sanitizers perfect for home, hotels, hospitals or other public places.
- The unit has 2 display LEDs to indicate ON/OFF status and dispensing in progress
- Ultrasonic sensor to ensure zero touch, completely disinfect both hands at once. Just put your hand under the sensor to start the soap dispenser, You can avoid cross infection without touching the soap dispenser, Standard 220 V AC from wall socket. Voltage actuation and surge protection is built in.

### SANITIZER STAND FOOT OPERATED



- Cost effective solution for Hand sanitization
- Manual machine with foot operated lever to dispense the sanitizer
- Easy to operate with any size of Dispenser Pump Bottle
- Best Suited for Offices, Clinics, Hospitals, Schools, Banks, Courier Services, etc.
- Hand Touch Free, Long Lasting, Good Quality, Very Sturdy and Maintenance Free
- Sanitizer bottle can be easily replaced by rotating knob
- Stand can be easily mounted at any surface within seconds
- Powder coated for long life and sturdy for maximum use at many places



## PRODUCT RANGE

### PPE KIT



Package contains

- Coverall Suit with Hood, Face mask, Hand gloves, Shoe leggings, Hood cap
  - Coverall suit 90 GSM, Sitra Certified for fabric and seams – penetration test
  - 3 ply masks or N 95 Masks as required
  - Shoe Covers Elongated till calf length providing 100% overlapping with coverall suit
  - Nitrile Gloves Better sensitivity and better elasticity
  - Chemical resistant Provide precise grip in both dry and wet conditions
  - Biohazard Bags 100% Recyclable Compatible with international standards
- PPE kit can be customized also in terms of sizes, more contents etc.

### SANITIZING LIQUID



- Sodium Hypochlorite Solution - SANITIZING Liquid
- Content 5% / 10 % of Sodium Hypochlorite.
- all-purpose cleaner for sanitation purpose
- House hold cleaner
- for sanitization use in office, home ,Hotel etc.
- germs remover
- Cost effective solution for sanitization
- Packing – available in 1 Liter / 3 liter/ 5 Ltr and 10 Ltr

### MANUAL SANITIZES SPRAYER



- Manual sanitizer sprayer machine of various capacities 4/6/8/10/12 and 16 Ltr
- Unbreakable Strong Base
- Capacity: many options available for bottle capacities, ISI Marked Sprayer
- Conventional Design and Strong Tank with brass Pressure Chamber
- Continuous Mist Spray and cozy and Adjustable Cotton Belt
- Strong Blow Molded Tank, effective solution for sanitization
- Manual machine with operated lever to dispense the sanitizer
- Easy to operate with any size of Dispenser Pump Bottle
- Best Suited for Offices, Clinics, Hospitals, Schools, Banks, Courier Services, etc.

### HANDRUB SANITIZER



- Ethanol-based Hand Rub Sanitizer and Disinfectant – various sizes
- 70 - 80% Ethanol based, Rub-in Hand Disinfectant and Sanitizer. Effective against 99.9% Germs, Viruses and Bacteria
- World Health Organization (WHO) recommended hand-rub formulation
- 0.2 Micron filtered Rub-in Hand Disinfectants & Sanitizers are pure, hygienic and effective during routine, high risk and outbreak situations
- It is ready to use, does not require any water, or towel drying, and can be used conveniently and as frequently as needed, Best Suited for Offices, Clinics, Hospitals, Schools, Banks, Courier Services, etc. Available in 10 ml, 50 ml, 100 ml, 500 ml, 1 liter and 5 Ltr etc.

## PRODUCT RANGE

### FOREHEAD THERMOMETER



- Digital LCD No Contact Infrared Thermometer-for Adult and kids
- High accuracy and fast- with the most advanced infrared technology, can acquire the aim temperature quickly and displaying the result within 0.5 second.
- Colour display- color backlight: green backlight for normal, yellow backlight for slight fever, red backlight for high fever
- Safety-it do not cause any un comfort and breed bacteria with the creative design non-contact any part of the body, especially for baby and elder, Easy to use- quiet mode tympanic scanning with the push of a button allows you to take instant scans of your baby without a beeping noise that may wake them
- Package: 1 \* Non-contact IR Thermometer 1 \* Manual (English) 2 \* Batteries.

### DISPOSABLE HEAD/SHOE COVER



- Disposable Stretchable Caps/Head Cover
- Best to Keep your hands and head clean while doing work. Best for: Home and workspace
- Simple to Use and Comfortable to wear with fresh new and high quality
- Material : Non-Woven fabric disposable Shoe Cover
- Usage : Hospital, clinic, food industrial, clean room, laboratory, home etc.
- Eco-Friendly, Biodegradation
- Pack of 25/50/100 Pcs
- High-quality elastic, durable double elastic band bundle
- tear-resistant, not easy to break More breathable and anti-skid

### WORKSPACE SAFETY SHIELD



- Especially designed to easily adapt existing workspaces for the post COVID-19 work environment.
- SAFETY SHIELD provides physical social distancing barriers for hospitals, pharmacies, grocery stores, check stands, and other point-of-purchase locations.
- Helps to prevent the spread of respiratory diseases while maintaining customer experience.
- Ideal for cashier, front desk, and hospitality employee additional protection.
- Food Service & Hospitality, Office / Reception Desks, Medical / Dental Offices, Booth Partitions

### PORTABLE SANITIZER



- Usable for Sanitizing Mobiles, Keys, Car, Wallet, Notes, Remotes, Watches , Laptops, Currency
- Features: The handy facial steamer can be used as a portable charger, power capacity is 150 mAh which can charge your cell phone for more than 1 hour.
- Portable size for everyday care: small size and light weight. You can conveniently put it in your backpack, sanitize it all day, This can also be used as a Water Mist to refresh the Skin
- Cooling method for sunburn: Producing ultra-fine and nano-cold fog will help soothe sunburn areas and make burns heal faster. In the hot summer, it helps to keep the skin hydrated, thus offsetting the dryness of sunburn.



## PRODUCT RANGE

### URRÀ KIT, IMMOBILIZATION KIT




urrà kit collects equipment for the treatment of fractures and sprains in a small bag with compartments made of strong material.

It contains three semi-rigid splints in different sizes, one instant ice pack, universal bandage scissors and self-adhesive flexible bandage.

It contains

- three Urrà Splints in different sizes
- one Spencer Ice Pack instant ice
- Universal bandage scissors
- white self-adhesive flexible bandage (10 cm x 6,3 m) [(3.93 x 248.03 in)]
- bag made of strong material, complete with compartments to contain all that is necessary

### HIRVIN - TRANSFER SHEET



- Scrupulous controls guarantee that every sheet responds to high standards of quality.
- In fact the Spencer logo is a worldwide symbol of high quality.
- during the construction of the Hirvin injury transportation sheets a lot of attention was given mainly to the ergonomics aspect, considering the actual structural
- dynamic conditions as during use by the operator. Invisible details that express its robustness and manageability.

### B-LIFE REUSABLE MANUAL RESUSCITATOR



- easy to remove pressure limiting valve
- transparent oxygen reserve valve
- disposable reservoir bag
- face mask
- PEEP valve
- Available in the following versions:
- adult: for patients over 30 kg [66 lbs]
- paediatric: for patients between 7 and 30 kg [between 15 and 66 lbs]
- Single packaging in cardboard box printed using a four-coloured process.

### TRANSPORT VENTILATOR, 190 NXT



- Dimensions: 330 x 177 x h210 mm, Weight: 3.65 kg
- Flow: from 2 to 20 L/min, Frequency: from 5 to 40 BpM
- I/E connection: 1:2
- Oxygen concentration: 60/100 %
- Battery duration: 3 h
- Alarms: high pressure, low pressure, apnea, pneumatic supply, low battery
- Alarms sound pressure level: 80 ± 10 db
- Audio mute time: 100 s
- Overpressure valve adjustment: da 20 a 60 cmH2O
- Manometer: from -20 to +80 mbar
- Pneumatic connection standard: UNI



Registered Office :

**Spencer India Technologies Pvt. Ltd**

Plot No. 477, Ground Floor, Sector 26, Bharthal, Dwarka, New Delhi 110077,

Contact : +91 99715 86305, [info@spencer-india.com](mailto:info@spencer-india.com)

[www.spencer-india.com](http://www.spencer-india.com)

