

The value of a choice is measured in the future.



Skid Series

Introducing the ultimate expression of Spencer's efficiency, design and innovation.

spencer.it

SPENCER[®]

SELF-BRAKING SLIDES
INTUITIVE USE
MAXIMUM COMFORT

SKID SERIES

Evacuation stair chair

Pro Skid evacuation chairs are designed with a caterpillar reinforced rubber belts, characterized by strong adherence to the surface together with easy maneuverability. The use of passive caterpillar belts allow an elevated level of security while they can effortlessly glide downstairs.

Unique features

- PVC sheet easily cleanable equipped with quick release buckles to adjust the tension

Pro Skid caterpillar belts can be easily assembled and replaced.



The seat, in PVC, is reinforced and adjustable, easy to clean and to replace in case of wear.

